



CONTROL PANELS

LGK-2 Digital



Microprocessor control panel with easy-to-read **digital display** showing inside menu, temperature, set point, date and hour, model of appliance.

Equipment settings and switching off only by means of **password**.

Continuous **temperature monitoring** with set point and alarms shown on display.

Visual and acoustic alarms through buzzer sound, general alarm warning light, specific alarm warning light and icon on display.

Alarms included: high and low temperature (adjustable threshold), open door, condenser efficiency, damaged probes, power failure, run-down battery.

Back up battery (optional) to show temperature on display and generate alarms in case of power failure.

Remote contact to connect the appliance to an external alarm.

GSM module (optional) to send SMS alarm messages showing date, time, cabinet temperature and type of alarm. Possibility to send SMS message to check correct functioning of cabinet.

LAN port (optional) to connect one or more appliances via ethernet on PC for monitoring correct functioning.

USB port to download temperature data and alarms or to reset the appliance.

Data logger to record temperatures of last 30 days. Download data on USB pen drive.

Memory for last 20 alarms shown on display.

